

ShowCase® Image Center

Study Storage, Local and Remote Image Viewing

STORAGE AND RAPID ACCESS TO 200,000+ IMAGING STUDIES

ShowCase Image Center runs on affordable hardware, from a laptop to a server to a VM, and will accept and store DICOM studies sent from your ultrasound machine or other imaging modality. Quickly locate studies using the powerful search features of the Image Center.

VIEW STUDIES USING SHOWCASE PREMIER OR VIEWER FOR IPAD

Access and view studies locally, in office, or remotely, via VPN, using the ShowCase Premier viewer or Viewer for iPad (ViP).


SHOWCASE PREMIER



ShowCase Image Center combined with ShowCase Premier allows you to view images in full quality DICOM on local and remote PCs.

ShowCase Premier runs on desktop or laptop PCs with Windows 7, 8.1, 10, 2008r2, 2012, or 2012r2. Recommended minimum of 4 GB RAM, 4 GB free disk space, high quality display.

SHOWCASE VIEWER FOR IPAD (ViP)



ShowCase ViP provides portable, non-diagnostic study viewing, locally or remotely. Ideal for showing images to patients or colleagues.

View single and multiple image layouts and play clips. Interact with displays using familiar iPad gestures.

ShowCase ViP is not FDA cleared for diagnostic use. ShowCase Viewer for iPad runs on iPad 2 or later and requires ShowCase Image Center software to convert and transfer images.

DESIGNED TO SCALE TO YOUR NEEDS

Small clinic looking for an affordable solution? No problem: install Image Center and Premier together on an affordable PC. Need to scale to a 20TB, 30 million image archive? Again, no problem: just provide the necessary hardware and disk space and you are good to go.

ALL DATA STORED IN NATIVE DICOM FORMAT

Image Center stores your data as it comes from your imaging machines, using the DICOM standard.

Caution: Federal law restricts the sale or use of ShowCase Premier by or on the order of a physician. ShowCase is a trademark registered in the U.S. Patent Office to Trillium Technology, Inc, Ann Arbor, Michigan.